



spa

via Croazia, 2 - 33100 UDINE - ITALY
tel. +39 0432.621621
fax. +39 0432.621620
www.chinesport.it - chinesport@chinesport.it

STAMPA LISTA RICAMBI / SPARE PARTS LIST

LT2B12.-001 - TILT TABLE 2BT

S-N.:

Listino:



Cliente: -

Vs.Rif.:

Ns.Rif.:

Manuale:

SECTION POSITION	INTERNAL CODE	DESCRIPTION	Foto/Image	Q.TY	Pris per enhet. SEK. Exkl. moms	Tipo ric./ Type part
-	85.24.01.0547-	CHINESPORT LABEL MM 92X42		PZ 2	30	-
-	85.24.01.0206-	COPRIDADO ESAG.13 - DADO M 8		PZ 4	5	-RI
-	85.24.01.0308-	HOLECOVER BLACK DIAM.14 X 8	image not available	PZ 4	5	-
-	85.24.01.0525-	TRIANGLE ADHESIVE LABEL FEET	image not available	PZ 3	10	-
-	85.24.01.0402-	HEX WRENCH MM 13/10		PZ 1	10	-



spa

via Croazia, 2 - 33100 UDINE - ITALY
tel.+39 0432.621621
fax.+39 0432.621620
www.chinesport.it - chinesport@chinesport.it

STAMPA LISTA RICAMBI / SPARE PARTS LIST

LT2B12.-001 - TILT TABLE 2BT

S-N.:

Listino:



Cliente: -

Vs.Rif.:

Ns.Rif.:

Manuale:

SECTION POSITION	INTERNAL CODE	DESCRIPTION	Foto/Image	Q.TY	Pris per enhet. SEK. Exkl. moms	Tipo ric./ Type part
-	85.24.01.0175-	HEX WRENCH MM 6	image not available	PZ 1	10	-RI
-	85.00.00.0017-000	TILT TABLE PVC LABEL WHITE STICKER		PZ 2	20	-
-	85.00.00.0028-	WARNING LABEL TILT TABLE MM 75X25 WHITE		PZ 2	20	-
P1	PI009.1A060-	PADDED HEADREST TILT WITH HOLE	image not available	PZ 1	2500	-
P2	PI010.1A060-001	PADDED CENTRAL TILT	image not available	PZ 1	1900	-



via Croazia, 2 - 33100 UDINE - ITALY
tel.+39 0432.621621
fax.+39 0432.621620
www.chinesport.it - chinesport@chinesport.it

STAMPA LISTA RICAMBI / SPARE PARTS LIST

LT2B12.-001 - TILT TABLE 2BT

S-N.:

Listino:



Cliente: -

Vs.Rif.:

Ns.Rif.:

Manuale:

SECTION POSITION	INTERNAL CODE	DESCRIPTION	Foto/Image	Q.TY	Pris per enhet. SEK. Exkl. moms	Tipo ric./ Type part
P3	PI012.1A028-	PADDED CALF LENGHT TILT	image not available	PZ 2	950	-
-	82.00.00.0037-	BATTERY 12V 2.7AH	image not available	PZ 2	750	A-R1
-	80.23.00.0134-	BASIC BRAKING WHEEL Ø125		PZ 1	250	-
-	80.23.00.0135-	BASIC STEERING WHEEL Ø125		PZ 1	250	-
-	80.23.03.0029-	BATTERY CHARGER		PZ 1	1250	-RI



via Croazia, 2 - 33100 UDINE - ITALY
tel.+39 0432.621621
fax.+39 0432.621620
www.chinesport.it - chinesport@chinesport.it

STAMPA LISTA RICAMBI / SPARE PARTS LIST

LT2B12.-001 - TILT TABLE 2BT

S-N.:

Listino:



Cliente: -

Vs.Rif.:

Ns.Rif.:

Manuale:

SECTION POSITION	INTERNAL CODE	DESCRIPTION	Foto/Image	Q.TY	Pris per enhet. SEK. Exkl. moms	Tipo ric./ Type part
-	80.23.03.0101-002	CBJ2 CONTROL BOX		PZ 1	1950	-RI
-	80.23.01.0428-	GAS SPRING TL=245,5MM C=80 800N		PZ 1	350	-
-	80.23.00.0004-	HEIGHT ADJUSTMENT ACTUATOR		PZ 1	1950	-
-	80.23.03.0431-	SLS NOR.APE. DINFEMM.100CM JACKMAS.130CM		PZ 2	750	-
-	01132-001	SPARE BATTERY		PZ 1	1950	-



via Croazia, 2 - 33100 UDINE - ITALY
tel.+39 0432.621621
fax.+39 0432.621620
www.chinesport.it - chinesport@chinesport.it


STAMPA LISTA RICAMBI / SPARE PARTS LIST

LT2B12.-001 - TILT TABLE 2BT

S-N.:
Listino:



Cliente: -
Vs.Rif.:
Ns.Rif.:
Manuale:

SECTION POSITION	INTERNAL CODE	DESCRIPTION	Foto/Image	Q.TY	Pris per enhet. SEK. Exkl. moms	Tipo ric./ Type part
-	80.23.00.0198-	WHEEL Ø100		PZ 2	250	-
-	85.24.01.0729-	PISTON GN 613-8-AK M16X1.5	image not available	PZ 2	175	-
70.00.00.0006/1	80.22.00.0544-001	BOARD RIGHT TILT TABLE STRATIF. WHITE	image not available	PZ 1	295	-
70.00.00.0006/2	80.22.00.0545-001	BOARD LEFT TILT TABLE STRATIF. WHITE	image not available	PZ 1	295	-
-	85.24.01.0555-	FERRULE 80X30 BLACK	image not available	PZ 2	20	-



via Croazia, 2 - 33100 UDINE - ITALY
 tel.+39 0432.621621
 fax.+39 0432.621620
 www.chinesport.it - chinesport@chinesport.it


STAMPA LISTA RICAMBI / SPARE PARTS LIST

LT2B12.-001 - TILT TABLE 2BT

S-N.:
 Listino:



Cliente: -
 Vs.Rif.:
 Ns.Rif.:
 Manuale:

SECTION POSITION	INTERNAL CODE	DESCRIPTION	Foto/Image	Q.TY	Pris per enhet. SEK. Exkl. moms	Tipo ric./ Type part
-	80.23.03.0430-	ACTUATOR FOR STATIC POSTURE 140KG		PZ 1	1950	-

NESSUNA IMMAGINE DI ESPLOSO ASSOCIATA / EXPLODED VIEW NOT AVAILABLE

FINE LISTA / END OF LIST