

Spare parts list MOTomed loop

Order no. 260.010/020/030/039/040

260.060/070/080/090

260.100/110/120

valid from January 2021 The prices indicated are subject to change!

Image No.	Name	Legend: A) Repairs can be carried out by customers, authorized partners and in the factory B) Repairs can be carried out by authorized partners and in the factory C) Repairs can only be carried out in the factory	Quantity Per Device	ID no./ Order no.	Spare part price [SEK / pc.] without VAT
Basic module (Sheet 1/13)					
		Execution	1		
1	Transport wheel right, hub closed, set	A	1	100020008MT	540
2	Transport wheel left, hub open, set	A	1	100020009MT	540
3	Power supply plug-in unit left, set	B	1	100020010MT	3 465
4	Plug-in unit with wheel right, set	A	1	100020011MT	625
5	Standing feet back, set	B	1	100020012MT	265
6	Standing feet front, set	A	1	100020013MT	135
7	Screws M3x10 IS7045 4.8 [galvanized], set	A	1	100020185MT	27
8,9	Covering caps front leg, set	A	1	100020014MT	150
Motor section (Sheet 2/13)					
			1		
1	Arm and leg trainer motor (GR63) for loop edition	B	1	300983101MT	4 675
2	Motor mounting, set	B	1	100020015MT	480
3	Power controller, set with software	B	1	100020016MT	5 380
4	Data bus line type F	B	1	300729101MT	430
5	Load resistor, set	B	1	100020017MT	150
6	Cable tie, set	B	1	100020018MT	45
Adjustment unit (Sheet 3/13)					
			1		
1	Bushing adjusting shaft [LO-VE]	B	1	100019365MT	180
2	Adjusting shaft [LO-VE]	B	1	100020108MT	1115
3	Guide tube adjusting shaft (until H04-LP-xxx) [LO-VE]	B	1	100014787MT	610
	Guide tube adjusting shaft (since H05-LP-xxx)	B	1	100019364MT	610
4	Twist protection adjusting shaft coated	B	1	100019360MT	855
5	Gear wheel adjusting shaft [black] [LO-VE]	B	1	100013271MT	39
6	Gear wheel servo motor [black] [LO-VE]	B	1	100013272MT	43
7	Servo motor 180°	B	1	300984102MT	2 200
8	Sensor, set	B	1	100020023MT	560
9	Line for position sensor servo	B	1	301011100MT	135
10	Supply of servo motor, set	B	1	100020025MT	605
11	Servo control line	B	1	301021100MT	135
12	Pan-head screw ISO 7045 M3x6, set	B	1	100020247MT	10
13	Pan-head screw ISO 7045 M3x10, set	B	1	100020248MT	10
14	Dowel pin POM ø4x10 [LO-VE], set	B	1	100020249MT	27
15	Disc ISO 7089 - 3,2, set	B	1	100020250MT	2
16	Rubber washer 8x3,2x1,0, set	B	1	100020251MT	10
17	Safety ring set, DIN 471 A10x1 (until SN H04...)	B	1	100020252MT	13
	Safety ring set, A8x0,8 (since SN H05...)	B	1	100020253MT	13
Arm trainer module (Sheet 4/13)					
			1		
1	Arm trainer axle, set	B	1	100020027MT	4 037
2	Arm trainer locking, set	B	1	100020028MT	515
3	Position sensor arm, set	B	1	100020029MT	565
4	Line position sensor arm trainer [LO-AM]	B	1	300982100MT	530

Spare parts list MOTomed loop

Order no. 260.010/.020/.030/.039/.040
260.060/.070/.080/.090
260.100/.110/.120

valid from January 2021 The prices indicated are subject to change!

Image No.	Name	Legend: A) Repairs can be carried out by customers, authorized partners and in the factory B) Repairs can be carried out by authorized partners and in the factory C) Repairs can only be carried out in the factory	Quantity Per Device	ID no./ Order no.	Spare part price [SEK / pc.] without VAT
Leg trainer axle (Sheet 5/13)					
		Execution	1		
1	Leg trainer axle, set	B	1	100020031MT	630
2	Line position sensor leg trainer [LO-BV]	B	1	300981100MT	545
Operating panel /transport handle (Sheet 6/13)					
			1		
1	Operating panel 7", set	B	1	100020033MT	555
2	Carrier card, set (with software)	B	1	100020034MT	6 355
3	FFC line	B	1	300986100MT	80
4;6;7	Swivelling set	B	1	100020037MT	270
5	Lower casing	A	1	100011538MT	135
6	Screws STS KN1032-new 30x16-H1-A2-1.4567, set	A	1	100020186MT	53
8	Lamella caps GL22x0.8-3 [black], set	A	1	100020187MT	42
9	Transport handle, set	B	1	100020041MT	705
Casing (Sheet 7/13)					
			1		
1	Plastic casing, set	A	1	100020042MT	1 980
2	Fixation, set	A	1	100020043MT	60
3	Clamping piece, set	B	1	100020044MT	25
Adjusting frame (Sheet 8/13)					
			1		
1	Locking pin D6 stainless steel M10x1-Type B	A	1	100010749MT	300
2	Control knob cone shaped	A	1	100014205MT	425
3	Handle mount, set	B	1	100020047MT	250
4	Position sensor board leg, set	B	1	100020048MT	500
5	Leg trainer mount coated	B	1	100012976MT	1 255
6,7	Screws, set	B	1	100020050MT	80
8	Hexagon head screw ISO 4017 - M8x70 - 8.8 -	B	1	100001453MT	19
SoftGrip handles (Sheet 9/13)					
			1		
1	Handles (pair), set	A	1	100020052MT	1 850
2	Clamping caps SW 13-1-G8 [black], set	A	1	100020188MT	76
3	Eccentric lever clamping system M5x20 adjustable	B	1	100012699MT	980
4	Plastic spacers ø25xø12x29, black, set	B	1	100020189MT	95
5	Pedal crank with parallel adjustment with arm	B	1	100015663MT	2 110
6	Safety ring DIN 471 A10x1, set	B	1	100013285MT	20
7	Handle mount, set	B	1	100020058MT	385
	Plastic-coated handles (pair) with quick release clip syst.	A	1	100020059MT	3 565
Plastic-coated foot shells (Sheet 10/13)					
			1	not for 260.100/.110/.120	
1	Plastic-coated safety foot shells flat complete	A	1	100017029MT	5 545
2	Clamping caps SW 13-1-G8 [black], set	A	1	100020188MT	76
3	Plastic plugs SFL 32x0,8-2,5 [black], set	A	1	100020190MT	57
SoftGrip foot shells (Sheet 10/13)					
			1	only for 260.100/.110/.120	
1	Foot shells (pair) flat complete mounted	A	1	100016911MT	4 820
2	Clamping caps SW 13-1-G8 [black], set	A	1	100020188MT	80
3	Plastic plugs SFL 32x0,8-2,5 [black], set	A	1	100020190MT	57

Spare parts list MOTomed loop

Order no. 260.010/.020/.030/.039/.040
260.060/.070/.080/.090
260.100/.110/.120

valid from January 2021 The prices indicated are subject to change!

Image No.	Name	Legend: A) Repairs can be carried out by customers, authorized partners and in the factory B) Repairs can be carried out by authorized partners and in the factory C) Repairs can only be carried out in the factory	Quantity Per Device	ID no./ Order no.	Spare part price [SEK / pc.] without VAT
Gas spring support (Sheet 11/13)					
			1	for 260.039	
1;3; 4;5	Fixation set	B	1	100020067MT	265
2	Gas spring GF6-ABAB30 / 31-200	B	1	100014260MT	2 970
Pedal radius quick adjustment / QuickFix (Sheet 12/13)					
			1		
1.1	Plastic-coated foot shells (pair), pedal rad. fixation straps	A	1	100017500MT	10 090
1.2	Plastic-coated foot shells (pair), pedal rad. disinf. fix. straps	A	1	100018024MT	10 978
2	QuickFix foot fastening system (pair)	A	1	265.004	5 450
3	Dipped caps 10x3x40mm [orange RAL 2004]	A	1	100020191MT	43
4.1	Leg guides with plastic-coated calf shells, fastening straps	A	1	303.000	1 950
4.2	Leg guides with plastic-coated calf shells, disinfectable straps	A	1	186.000	2 500
5	Wing screws M4x8, set [black]	A	1	100020192MT	46
Transmission axle (Sheet 13/13)					
			1		
1	Transmission axle complete, set	C	1	100020197MT	5 745